# NUVOLINA LED LIGHT

MANUALE D'USO

USER'S GUIDE


G.COMM

**INDICE**

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**CONTENTS**

SYMBOLS USED IN THE MANUAL AND IN THE PRODUCT..... 7  
 SAFETY WARNINGS..... 7  
 INTRODUCTION..... 7  
 TECHNICAL DESCRIPTION..... 8  
 PLACEMENT..... 9  
 INSTALLATION..... 9  
 ELECTRICAL WIRING..... 10  
 USE AND MAINTENANCE..... 10  
 EMERGENCY KIT..... 10  
 STORAGE CONDITIONS..... 10  
 REPAIRS..... 10  
 WARRANTY..... 12  
 WARNING RAAE..... 12

**SYMBOLS USED IN THE MANUAL AND IN THE PRODUCT**

- Follow instructions for use
- Prohibition (this icon identifies what to prevent doing not to cause damage to the device)
- Warning (this symbol warns to be very careful in order to avoid any situation that could cause damage to the device)
- RAEE symbol
- Ground
- Alternating current

**SAFETY WARNINGS**

**Before using the operating light system, please read the instructions contained in this manual and please pay attention to the following indications:**

- the ceiling has the intended use of workplace lighting and should only be used for this purpose;
- the installation of the device should only be performed by qualified personnel;
- the lamp must be installed and powered according to the power requirements present in this manual and with electrical equipment which complies with the standards applicable to rooms used for medical purposes;
- do not perform any maintenance or cleaning on the lamp when the power is inserted;
- do not insert any objects into the slots of the ceiling light that may touch dangerous points under voltage;
- ensure cleaning and disinfection of ceiling light's surface as specified in this manual;

- do not cover the device in order to ensure ventilation;
- use of vapor or spray liquids near of the lamp prevent them from entering the housing of the light fixture itself.

**INTRODUCTION**

**Thank you for purchasing a ceiling light of "NUVOLINA Light led" series.**

Designed according to three main principles: aesthetics, functionality and simplicity, with particular attention to design and materials. The result is a new ceiling light, with smooth and simple lines, extremely functional and solid. "NUVOLINA Light led" provides superior optical performance level to meet every professional dentist's need.

The series of professional ceiling lights "NUVOLINA Light led" is made with materials that permit easy cleaning and prevent bacterial growth. The elegant design and structure without openings make it perfect for use in medical and maintenance free environments.

**FUNCTIONS**

Thanks to the use of the most innovative LED light source this series of ceiling lights ensures ideal lighting levels, in compliance with the UNI EN 12464-1, in three areas of the treatment room:

- at least 1500 lux on the patient's face
- at least 1000 lux in the work area
- 500 lux in the transition zone.

• **High color rendering index (> 90):**  
 G.Comm's ceiling light "NUVOLINA" series ensure an high color rendering index and provides excellent color vision in the mouth.

• **Light sources of latest generation:**  
 equipped with LEDs that can provide a colour temperature (K):

- 4500 or 5500 K for NUVOLINA Light led
- from 3000 to 5700 K for NUVOLINA Light led Tunable White.

• **Long light lifetime :**  
 the LED technology allows a long lifetime of the light (around 50.000 hours)

**For any need, contact our Customer Service +39 039 6060420**

**G.Comm disclaims any liability arising from in**

**SAFETY WARNINGS**

**Before using the operating light system, please read the instructions contained in this manual and please pay attention to the following indications:**

- the ceiling has the intended use of workplace light-

- ing and should only be used for this purpose;
- the installation of the device should only be performed by qualified personnel;
- the lamp must be installed and powered according to the power requirements present in this manual and with electrical equipment which complies with

**TECHNICAL DESCRIPTION**


Figure 1

	<i>Nuvolina Light led Single Color</i>	<i>Nuvolina LED Tunable White</i>
Ceiling height [m]	2.4 - 3	2.4 - 3
Energy consumption [W]	300	300
Maximum luminance [lux] at 850 mm and 1400mm	7000 - 3500	7000 - 3500
Luminous flux [lm]	20700	20700
Colour temperature [K]	5700	from 3000 to 5700
Color rendering index [%]	>90	>90
Dimensions [mm]	1176 x 510 x 86	1176 x 510 x 86
Supply voltage	220/240 Vac 50/60 Hz	220/240 Vac 50/60 Hz
Weight [Kg]	15	15

PLACEMENT


Figure 2

It's suggested to place the center of the lamp at the workplace's center as shown in fig. 3.

INSTALLATION

**! WARNING**

The device must only be installed by qualified technicians.

The installation and all maintenance operations must always be carried out without power supply. Before installing the device control that electrical data of your system are in accordance with those in this manual.

- Secure the template with holes position (D) with the adhesive tape to the ceiling;
- By the means of a drill (A) carry out the holes (tip diameter 12 mm) as indicated on the template set previously;
- Install the conical rosettes (B) using an appropriate fastening system (i.e. butterfly blocks (C), expansion or chemical plugs) depending on the ceiling type;
- Insert the cable seal (F) in the final (E) and screw it clockwise into the conical rosettes (B);
- Place the lamp on the ground still inside the packaging in correspondence with the cables (F) (cable length is 3 m);
- Insert the cable into the terminals (H), using the following procedure: unscrew the outer ring and push the clamp down sliding in the cable;
- Begin to secure the cable on the short side of the lamp;
- Raise the lamp using F1 side as starting point;
- Lift the lamp on the side leaning on the ground and secure the cable to the terminals as explained in the previous step;
- Lead the lamp to the desired height, gradually pulling the cables of a short side and then the others;
- Once the desired height is reached, place a level on the upper side of the lamp and adjust the strain of the cables;
- Cut the remaining cable paying attention to the

electrical connection of the lamp holder. Finally close the two sides of the ceiling light with the two chassis provided.


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7


Figure 8


Figure 9

ELECTRICAL WIRING

The electrical connection of the ceiling light is done by connecting the 3 - pin cable as follows:

- the transparent cable identified by the blue wire must be connected to the neutral of the interrupted branch of the power supply line;
- the transparent cable identified by the brown wire must be connected to the phase of the interrupted branch of the power supply line;
- the transparent cable identified by the green wire should be connected to the ground of the power supply line;

USE AND MAINTENANCE

The ceiling lights "NUVOLINA Light led" have been designed with particular attention to hygiene to allow for easy disinfection.

It will be very simple to clean your lamp following these brief tips:

- **Light diffuser:**  
Use a soft cloth moistened with water caressing the surface, and dry with soft cloth, better if antistatic.

**! AVOID TO CLEAN UP DIFFUSERS WITH ALCOHOL OR ABRASIVE SUBSTANCES.**

- **Lamp framework**  
For the usual cleaning, use a soft microfiber cloth dampened with water and dry with a soft cloth. For NUVOLINA LIGHT add a common household liquid detergent to the water and for persistent stains use an universal degreaser.

EMERGENCY KIT (optional)

The emergency kit, during electricity interruptions, generate a luminous flux in compliance with the current UNI EN 1838 standard, properly called: NO-panic light.

Part of the safety lighting of a building is intended to avoid panic and to provide light for identify safe escape routes.

After installation, the battery will reach full charge after 24hours.

The autonomy of the emergency luminous flux is more than 60 minutes every 12 hours of recharging according to regulations.

The Kit is equipped with a Green LED. When the green light is on and continues, the battery is charged and ready to work.

**! ATTENTION:**

The NO-panic light in the Nuvolina Light led does NOT replace, but integrates, the emergency lighting which, by current legislation, is necessary for the performance of the medical activity

of the dental office.

BATTERY REPLACEMENT

If the battery is no longer able to provide enough energy to power the emergency kit, replace the battery following the procedure:

**! ATTENTION:**  
Installation and all maintenance operations must always take place with the power supply disconnected.

1. Remove the side plastic by removing the 4 fixing screws;
2. Lift the central cover;
3. Disconnect the battery connector from the emergency kit by lifting the retaining clip.
4. Connect the new battery to the emergency kit using the appropriate connector;
5. Fasten the battery to the body of the lamp using a suitable clamp.
6. Replace the cover and the two side covers.
7. Tighten the 4 fixing screws.

**! ATTENTION:**  
Use only original battery packs supplied by G.COMM.

The battery is a high temperature Ni-Cd battery specific for emergency kits.

It supports a permanent charge for at least 4 years in high temperature environments (up to + 50 °).

Lifecycle over 500 recharges.

STORAGE CONDITION

Store in a closed, dry and protected place. Do not expose to heat stress and do not expose out of the following limits:

- Temperature: -10 ° C to 50 ° C
- Relative humidity: 10% to 90%
- Atmospheric pressure: 700 to 1060 hPa.

REPAIRS

G.Comm is never responsible for repairs carried out by third parties or unauthorized persons. Contact us directly to get the addresses of distributors and authorized technicians.

G.Comm is responsible, for the purposes of safety, reliability and performance of the equipment only if:

- Assembly operations, extensions, adjustments, modifications or repairs are carried out by authorized personnel;
- The electrical installation of the room complies with the requirements;

- The device is used in accordance with instructions provided. The product is covered by warranty for a period of **24 months** from the date of delivery to the end user.

#### WARRANTY

Warranty service must be performed by G.Comm or via an authorized repairer. Costs and risks of transport are borne by the customer. Warranty repair will be made only if the item to be repaired will be accompanied by invoice.

The warranty covers failure due to bad quality of material or manufacturing defects, in case of a claim founded the warranty will repair or replace free of charge.

There is no possibility of obtaining compensation for damages and / or interests. The guarantee is not valid, at the discretion of G.COMM, in case of tampering, damage, improper use, lack of maintenance or normal wear.

The neon is not covered by warranty.

#### RAEE WARNING

##### Indications about disposal

The European directives 2002/95/CE, 2002/96/CE and 2003/108/CE foresee the dental operating light in the category of electrical and electronic products that should be disposed separately from the municipal waste stream (RAEE).

The disused product shall be disposed separately to optimise the recovery and re-cycle of materials that compose the unit, in order to obtain a significant energy saving and to prevent potential negative consequences for environment and human health.

This crossed-out bin symbol is shown to all products that request to be disposed via designated collection facilities.

The user can hand over the disused dental operating light to his/ her city office or to the designated collection facilities appointed by his/her government or local authorities.

The illegal or inadequate disposal of disused product by the user implies the application of economic and/or administrative sanctions accordingly to the Law.

We reserve the rights to modify the specifications without prior notice.